

For GST-tag purification

Cell Lysis

- Lysis Buffer
- 40mM HEPES, pH7.5
- 0.5 M NaCl
- 10% glycerol
- 0.1% NP40
- Proteinase Inhibitor

For 10cm plate:

1. Remove medium from the plate.
2. Wash cells with cold PBS. (10ml)
3. Add 5ml lysis buffer to the plate. (or 2-5ml → increase protein conc.)
4. Leave on ice for 10min.
5. Shake the plate and transfer the cell lysate to a 15ml conical tube.
6. Centrifuge the lysate for 3-5min (2000rpm).
7. Recover the supernatant and store it at -80°C .

GST fusion protein purification:

1. Add 7.5 mg glutathione-agarose for each plate. Rotate at 4°C 60-90 min.
2. Load the mix into a mini-column. Wash with 5X5ml of buffer G plus 0.1% TritonX100 followed by 5X5ml of Buffer G.
3. Elute Gst fusion protein with Buffer G plus 10mM glutathione (0.031g/10ml), pH8.0 (glutathione is acidic so check pH of the final solution—need 12 μl 10M NaOH to bring 10ml of glutathione solution back to 8.0) in the following manner: resuspend beads in 1ml glutathione solution, incubate at room temperature for 10min, and then elute. Repeat elution 4 more times.

Buffer G

- 50mM Tris (pH8)
- 100mM NaCl
- 0.5mMEDTA
- 1mMDTT
- 10% glycerol

300mg
for
1 ml
lysates
out 5ml
total lysate
(1-2 150mm)
dish
293T