

PMAF Prog E 252-254

Induce 0.5d @ 30°C w/ 2µM IPTG 4hr.

Spin culture 6000 RPM 10'. Usually, 2 x 250ml.

pour off medium.

Add 5ml 1X lysis buffer to each bottle and resuspend pellets.

Freeze on dry ice.

Thaw @ 55°C.

Transfer to 50ml beaker

Rinse bottles with 1X lysis buffer and combine with previous 10 ml.

Final volume ~ 20ml.

Add Solid NaCl to 0.5M and DTT to 5mM final concentrations. 50% Duty cycle, Power @ 6

Sonicate with large tip 3 or 4 x 1.5 minutes.

Spin in ultra @ 23500 RPM (~75Kg) for 20 minutes.

Run gel to see if soluble, or not.

If Soluble:

purify on amylose resin and couple to Affi-gel 10

If Insoluble:

My insoluble fraction is usually in the same volume PBS as is the soluble fraction.

Spin PBS insoluble to re-pellet.

Save PBS soluble fraction if desired.

Assuming 0.5ml pellet volume, add 0.5ml 2x SDS sample buffer minus pms and CB.

Add another 4ml of 1X SDS sample buffer for a total vol. of 5ml.

(Everything should go into solution) Sonicate & boil if necessary

Dialyze 5ml against your coupling buffer, preferably o/v @ 4°C. Change buffer: RT w/ SDS

Couple to Affi-Gel 10.

Affi-Gel 10 Coupling:

Ideally, use 20-30 mg protein in a max. vol. of 4.5ml buffer for each ml of beads. \*Use MOPS, No Phos., HEPES, etc. NO TRIS buffers.

For my PMAF Prog E 252-254 I had ~20mg total protein solubilized in 5ml SDS sample buffer.

Take 1ml vol. of Affi-gel 10 in a 15ml screw-cap tube.

Wash 4x with 10 volumes cold dd H<sub>2</sub>O. Spin @ 900 RPM / 5'

Add 4ml/5ml soluble protein solution to beads.

Rock for 4 hours @ R.T. (Normally, at 4°C, but this will ↓ SDS)

After 4 hours, spin as above.

Save depleted protein solution. Run same on a gel.

Add 200µl 1M ethanolamine, pH 8 to bead pellet to block any active esters.

Rock @ RT 1hr (Again, usually @ 4°C).

After 1 hour, use TBS to transfer beads to a column, and to remove ethanolamine.

Beads can be stored in TBS + NaAzide @ 4°C.

# To Purify Antibody:

sera should be passed over MBP columns 2X for 2hr. before passing over protein affinity column. Procedure for all columns is the same.

Take column from cold.

Wash with 2 volumes TBS to remove azide.

Wash column with 25 volumes elution buffer to remove unbound protein. Do this before each use. For the first time after making the column, monitor the flow thru for protein content.

Wash column with ~5 volumes of TBS to remove acidic elution buffer.

Use TBS to transfer beads back into a screw-cap tube if the serum volume is too great for the column itself.

Spin @ 900 RPM / 5 minutes and remove TBS.

Add serum to beads.

Rock @ 4°C O/N for protein affinity column, or 2hr for MBP col.

Next day, transfer beads back to column using the serum itself to wash the tube.

Collect and save depleted serum for western.

Wash column with 25 vol. TBS to remove excess serum.

Column is ready to elute.

Eluate with fractions equal to bed vol.

Use 15-20 eppendorf's in a rack.

Add 200 $\mu$ l 1M Tris pH 9 to each tube.

Add 1ml elution buffer to top of column and collect. Repeat 15-20x.

Mix tubes well.

Eluted antibody must be neutralized immediately.

Wash col. with ~5 vol. PBS to remove acidic elution buffer (check with pH paper).

Read OD<sub>280</sub>.

100 $\mu$ l good fractions.

Microcon 30 to <0.5ml

Add original combined volume of Ca<sup>++</sup>, Mg<sup>++</sup> free PBS to concentrated antibody.

Microcon 30 to <0.5ml.

Add glycerol to 50% and freeze @ -20°C

## Example of dialysis/coupling buffer

50-100mM HEPES pH 8

0.1-0.25M NaCl

+ protease inhibitors.

2X  $\text{pH}$  buffer.

## Acidic Elution Buffer

0.5% Acetic Acid

0.15M NaCl

I usually don't add protease inhib.

	1x conc.	Stock	10ml	30ml	50ml
Naphos.	50mM	0.5M	2ml	6ml	10ml
MgCl <sub>2</sub>	5mM	1M	100 $\mu$ l	300 $\mu$ l	500 $\mu$ l
NaCl	0.2M	5M	0.8ml	2.4ml	4ml
BSA	0.1mg/ml	solid	2mg	6mg	10mg

$0.75\text{mg/ml} = 1 \text{ OD}_{280}$

$\rightarrow 50\text{ml}$

in the final  
+ protease inhibitor